direct emission lasers

# Superbium

LASER PLATFORM

ALL the Aesthetic Medicine in  
**ONE DEVICE**

 bios **LASER LINE**

# Company Profile

**BIOS** is an Italian Company that designs, manufactures and markets through a global distributors network, medical lasers and medical-aesthetic devices.

Founded in 1982, is **ISO 9001:2008, ISO 13485:2012 and 21CFR820 (USA)** certified. **BIOS** is grown up becoming a manufacturer company of hi-tech devices in compliance with government agencies requirements such as **FDA, ANVISA, GOST**. The technologies are designed to meet physician's needs and based on the worldwide daily experience and feedback of many experts and key opinion leaders, in order to verify the clinical validity of the methodologies and release the protocols ensuring the reproducibility and efficacy of the results.

The headquarters is located in **Milan** where workshops and training courses take place frequently. The medical devices market, strongly regulated and challenging, the **wide range of products and countries supplied** make the medical business very complex, where lean thinking, research of quality, innovation and voice of the customer are our key drivers to create value for our stakeholders.

**Strategic positioning** is oriented to three main business Areas.

**1>Medical Lasers**, for dermatology, surgery, dentistry where stands out:

**SUPERBIUM®**, the first and unique laser platform worldwide suitable for dentistry and aesthetic medicine through the combination of 7 lasers built-in the handpiece. **Superbium** is suitable for vascular lesions, wrinkles and scars reduction, hair removal, tattoo removal, skin photo-rejuvenation. At the same time can intercept the growing demand of dentists approached to the aesthetic medicine, supplying the useful tool of a laser to perform hard and soft tissues surgery.

**BIOSALEX PLUS YAG®** is an Alexandrite blended with Nd:YAG laser, the state of the art for hair removal; Is the most popular technology for hair removal, mostly in Middle East and Turkey, wherever in presence of dark skin type population. It has a long term and definitive results on unwanted hair and vascular lesions nevertheless the skin respect.

Complete the range of products the followings: **BIOPIXEL®** fractional and surgical CO<sub>2</sub>; **BIODEEPLUX®** diode lasers for ENT, EVLT, laser lipolysis.

**2>Medical-aesthetic devices** for cosmetic medicine and aesthetic applications on body and face such as: **BIOREVITAL RF MED®**, the 250W temperature controlled radiofrequency; **BIO SHAPE MED®** a multifunctional device for face skin tightening and body contouring based on bipolar RF, adipo-massage and IR.

**3>Physiotherapy devices: BIOREVITAL RF MED® Tekar Version**, useful for pain therapy, rehabilitation, reduce injury downtime.

Thanks to its steady and sustainable growth over the years, **BIOS** is nowadays positioned among the most reputable companies worldwide in terms of performances, reliability and quality of the devices. This is a basic skill to anticipate market trends with innovative products.

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# Superbium

## LASER PLATFORM

The new edge in dermatology and aesthetic medicine  
**ALL IN ONE LASER!**

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## Overview

### Panoramica

Starting today, performing wide range of treatments in the Doctor's clinic is no longer a huge investment in several different lasers. The solution is the innovative laser platform **SUPERBIUM**, an absolutely technical hedge that offers the possibility to perform dental surgery on both hard and soft tissues as well as cosmetic and dermal laser treatments. **SUPERBIUM** can work with 7 direct emission laser handpieces: **dermal Er:YAG 2940nm**, **dental Er:YAG 2940nm**, **fractional Er:Glass 1540nm**, **SP & LP Nd:YAG 1064nm**, **Q-switched Nd:YAG 1064+532nm**, **IPL** and **IR**.

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## Plug&Play

Its key to **versatility** lays on the **technological innovation** while separating the device management units (both electronic and hydraulic) from their optical unit where the laser is actually generated and the re-location of the latest inside the hand piece. This is why we don't speak anymore of "hand piece" alone but of "**laser hand piece**" each one of them has been equipped with a **PLUG&PLAY** connector which allows the Md to shift quickly between different sources available for use.

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## Assisted Mode

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Having experience in all laser applications is not easy and **SUPERBIUM** – with its countless disciplines where it can be used- actually puts on a hard test also for the most experienced MD. For this reason **SUPERBIUM** is gifted with the **assisted mode** which guides the MD through a series of questions and images towards the choice of laser settings more suitable for the treatment to perform according to the patient's photo type. The scientific sources collaborating with Bios **have studied more than 200 protocols** which have been incorporated within its software and are **now available for the operator** to use while working under its intuitive and **simple graphic interface**.

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
Save User's parameters  
Salvataggio parametri personali

## Adaptive electronic *Electronica adattiva*

**SUPERBIUM's real innovation stands in the advanced electronic technology.** Each laser source requires a specific power supply and custom fit hardware to achieve the maximum efficiency. That's why till today the common laser platforms on the market were powerless than stand alone lasers. **SUPERBIUM** breaks this trade off thanks to the innovative electronic technology **combines 3 different electronic hearts and power supplies** specifically developed for each laser source requirements: high frequency as Er:YAG needs, high energy as Nd:YAG and short pulse less than 12ns for Q-Switching. The result is the **amazing performances** of the laser handpieces comparable or **superior** than stand alone lasers with 2,5J for Er:YAG, 45J for Nd:YAG and 1800mJ for the Q-Switched.

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## Square Pulse *Impulso quadro*

**Advanced electronic allows to achieve square pulses.** Thus the laser energy delivered over the tissues remains **constant** during the entire pulse term, **keeping** the skin out of burnouts or side effects. Technology places **SUPERBIUM** at the edge of **reliability and safety**.

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


## Remote Service (TeamViewer®)

**Amazing remote service system:** SUPERBIUM display is not a simple touch screen, is a powerful computer windows based with all the potential of a such operative system. The SUPERBIUM control software is integrated with the popular **TeamViewer®** remote desktop software that allows the distributor and Bios service department to promptly perform the product assistance, remotely from the desk or the smartphone.

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TeamViewer® is copyright and property of TeamViewer GmbH

**Bios S.r.l.**  
 Via Guido Rossa 10/12 - 20090 Vimodrone (MI) - Italy  
 T +39 02 27304275 - F +39 02 27304276  
 info@biosgroup.eu - www.biosgroup.eu

MADE IN ITALY 0476

## Optional Accessories

Accessori optional

**Cart**  
 Carrello

**Trolley**

**Biocooling Air**

**Bag**

# Superbium

## **Nd** YAG SP&LP

45J x 1Hz (long pulse) - 1000mJ x 10Hz (short pulse) @ 1064nm

LASER LINE

Technology  
*Tecnologia*

SP	<b>1064</b> nm Wavelength	 Sapphire cooled tips	LP	<b>45</b> J Energy	<b>1</b> Hz Pulse Repetition Rate	<b>25</b> ms Pulse Width
			SP	<b>1000</b> mJ Energy	<b>10</b> Hz Pulse Repetition Rate	<b>200</b> µs Pulse Width

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


• Hair removal  
*Epilazione*

# The Multiplex Complex

Interaction between laser and blood vessels  
*Interazione laser vasi sanguigni*

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## Standard Accessories

*Accessori standard*

Ø3mm cooled sapphire tip  
 Ø6mm cooled sapphire tip  
 Ø9mm cooled sapphire tip

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Bios S.r.l. reserves the right to make changes on the technical product data.

This document does not represent an enforceable legal contract.

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
*Accessori optional*

- B** Ø1÷3.5mm tip

- C** 7x7mm 49 dots/cm<sup>2</sup> FRACTIONAL tip  
9x9mm 81 dots/cm<sup>2</sup> FRACTIONAL tip

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

### PEELING TREATMENT WITH QS LASER ASSOCIATED TO CARBON BASED LOTION TRATTAMENTO PEELING CON LASER QS ASSOCIATA A LOZIONE A BASE DI CARBONE

The 1064nm Q-switched laser, in synergy with a coal-based lotion; can perform a controlled and non-ablative peeling on the dermis surface and its remodeling. The laser beam heats the thin layer of skin fibroblast and collagen, thus providing a younger and healthier skin look. The effect is the skin tone improvement, narrow wrinkles, enlarged pores and reduction of epidermal or mixed melasma.

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



- Hair removal  
*Epilazione*

- Hair removal  
*Epilazione*

- Benign pigmented lesions  
*Lesioni pigmentate benigne*

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Operating view  
*Videata operativa*

Intended uses  
*Applicazioni*

- Skin rejuvenation through the soft tissues coagulation;
- Long pulse width for a deep penetration highly effective and powerful cooling;
- Advanced Contact Cooling™ technology to protect epidermis during the treatment;
- Skin tightening;
- Best results in combination with fractional non-ablative Er:Glass laser 1540nm, for a global face and body rejuvenation.

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# NIR: Near InfraRed technology

The NIR module is intended to emit energy in the infrared spectrum to provide topical heating for the purpose of:

- Elevating the tissue temperature for the temporary relief of minor muscle pain and joint pain and stiffness
- The temporary relief of minor joint pain associated with arthritis
- The temporary increase in local circulation where applied, and the relaxation of muscles
- May also help muscle spasms, minor sprains and strains, and minor muscular back pain.

## Skin Tightening and body Contouring using Near InfraRed technology.

The NIR successfully treats a wide range of areas, including:

- Skin tightening of body, face, and neck;
- Laxity of post-pregnancy abdominal skin;
- Laxity of buttocks and knees;
- Improving skin quality of décolleté.

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